





网版印刷电极



光滑碳电极




网版印刷电极




- 可抛弃式
- 可重复使用
- 可前处理后使用
- 不须研磨
- 可任意修饰
- 可大量使用
- 低价位
- 三电极(W.C.R)整合




Screen-Printed Electrodes 网版印刷电极

In-Built Three-Electrode Configuration

NO.#	TE100
Working Diameter/ Area	3 mm / 0.071 cm ²
Picture	



Single Working Electrode Configuration

NO.#	SE100	SE101	SE102
Working Diameter/ Area	5 mm / 0.196 cm ²	3 mm / 0.071 cm ²	1.5 mm / 0.018 cm ²
Picture			



CH Instruments

<http://www.chinstruments.com/>